

Product Datasheet

Human Phospholipase D2 ELISA Kit (orb1210195)

Description	Human Phospholipase D2, PLD2 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	PLD2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	105,987
Uniprot ID	O14939
Solubility (25°C)	160ng/ml
Sensitivity	0.29ng/ml
Entrez	5338
Expiration Date	Please enquire.