

Product Datasheet

Human Plexin C1 ELISA Kit (orb1210197)

Description	Human Plexin C1, PLXNC1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	PLXNC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	175,742
Uniprot ID	O60486
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	10154
Expiration Date	Please enquire.