

Product Datasheet

Human Plexin B1 ELISA Kit (orb1210202)

Description	Human Plexin B1, PLXNB1 ELISA Kit
Reactivity	Human
Range	50-10000ng/L
Tested Applications	ELISA
Target	PLXNB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	232,298
Uniprot ID	043157
Solubility (25°C)	12800ng/L
Sensitivity	25.84ng/L
Entrez	5364
Expiration Date	Please enquire.