

Product Datasheet

Mouse Vinculin ELISA Kit (orb1210249)

Description	Mouse Vinculin, VCL ELISA Kit
Reactivity	Mouse
Range	20-3500ng/L
Tested Applications	ELISA
Target	VCL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	116,717
Uniprot ID	Q64727
Solubility (25°C)	4000ng/L
Sensitivity	9.63ng/L
Entrez	22330
Expiration Date	Please enquire.