

Product Datasheet

Mouse Tumor Necrosis Factor-inducible Gene 6 Protein ELISA Kit (orb1210273)

Description	Mouse Tumor Necrosis Factor-inducible Gene 6 Protein, TNFAIP6 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	TNFAIP6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,924
Uniprot ID	008859
Solubility (25°C)	9600ng/L
Sensitivity	14.57ng/L
Entrez	21930
Expiration Date	Please enquire.