

Product Datasheet

Mouse Tumor Necrosis Factor Ligand Superfamily Member 14 ELISA Kit (orb1210280)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily Member 14, TNFSF14 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	TNFSF14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,338
Uniprot ID	Q9QYH9
Solubility (25°C)	3200ng/L
Sensitivity	7.65ng/L
Entrez	50930
Expiration Date	Please enquire.