

Product Datasheet

Mouse Thioredoxin-interacting Protein ELISA Kit (orb1210295)

Description	Mouse Thioredoxin-interacting Protein, TXNIP ELISA Kit
Reactivity	Mouse
Range	20-3500ng/L
Tested Applications	ELISA
Target	TXNIP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,363
Uniprot ID	Q8BG60
Solubility (25°C)	4000ng/L
Sensitivity	9.98ng/L
Entrez	56338
Expiration Date	Please enquire.