

Product Datasheet

Human Protein Tyrosine Phosphatase Type 4 A 3 ELISA Kit (orb1210338)

Description	Human Protein Tyrosine Phosphatase Type 4 A 3, PTP4A3 ELISA Kit
Reactivity	Human
Range	20-7000ng/L
Tested Applications	ELISA
Target	PTP4A3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,535
Uniprot ID	075365
Solubility (25°C)	8000ng/L
Sensitivity	11.98ng/L
Entrez	11156
Expiration Date	Please enquire.