

Product Datasheet

Rat Selenoprotein-P ELISA Kit (orb1210349)

Description	Rat Selenoprotein-P, SE-P ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	Se-P
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,083
Uniprot ID	P25236
Solubility (25°C)	3200ng/L
Sensitivity	8.13ng/L
Entrez	29360
Expiration Date	Please enquire.