

Product Datasheet

Mouse Serpin B8 ELISA Kit (orb1210350)

Description	Mouse Serpin B8, SERPINB8 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	SERPINB8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,150
Uniprot ID	O08800
Solubility (25°C)	4800ng/L
Sensitivity	9.11ng/L
Entrez	20725
Expiration Date	Please enquire.