

Product Datasheet

Rat Selenoprotein-S ELISA Kit (orb1210351)

Description	Rat Selenoprotein-S, SE-S ELISA Kit
Reactivity	Rat
Range	20-4500ng/L
Tested Applications	ELISA
Target	Se-S
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,385
Uniprot ID	Q8VHV8
Solubility (25°C)	4800ng/L
Sensitivity	10.13ng/L
Entrez	286900
Expiration Date	Please enquire.