

## Product Datasheet

# Mouse Serine Threonine-protein Kinase Tnni3K ELISA Kit (orb1210353)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Mouse Serine Threonine-protein Kinase Tnni3K, TNNI3K ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                          |
| <b>Range</b>               | 10-2000ng/L                                                    |
| <b>Tested Applications</b> | ELISA                                                          |
| <b>Target</b>              | TNNI3K                                                         |
| <b>Storage</b>             | 2-8°C                                                          |
| <b>Note</b>                | For research use only                                          |
| <b>Application notes</b>   | Sandwich                                                       |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                         |
| <b>MW</b>                  | 92,575                                                         |
| <b>Uniprot ID</b>          | <a href="#">Q5GIG6</a>                                         |
| <b>Solubility (25°C)</b>   | 2400ng/L                                                       |
| <b>Sensitivity</b>         | 5.34ng/L                                                       |
| <b>Entrez</b>              | <a href="#">435766</a>                                         |
| <b>Expiration Date</b>     | Please enquire.                                                |