

Product Datasheet

Rat Neurogranin ELISA Kit (orb1210354)

Description	Rat Neurogranin, NRGN ELISA Kit
Reactivity	Rat
Range	1-400ng/ml
Tested Applications	ELISA
Target	NRGN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,496
Uniprot ID	Q04940
Solubility (25°C)	480ng/ml
Sensitivity	0.63ng/ml
Entrez	64356
Expiration Date	Please enquire.