

Product Datasheet

Mouse Serine Protease Inhibitor Kazal-type 7 ELISA Kit (orb1210355)

Description	Mouse Serine Protease Inhibitor Kazal-type 7, SPINK7 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/mL
Tested Applications	ELISA
Target	SPINK7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,282
Uniprot ID	Q6IE32
Solubility (25°C)	32ng/mL
Sensitivity	0.023ng/mL
Entrez	408198
Expiration Date	Please enquire.