

Product Datasheet

Rat Haptoglobin ELISA Kit (orb1210356)

Description	Rat Haptoglobin, HPT ELISA Kit
Reactivity	Rat
Range	5-600ug/ml
Tested Applications	ELISA
Target	HP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,563
Uniprot ID	P06866
Solubility (25°C)	640µg/ml
Sensitivity	0.22µg/ml
Entrez	24464
Expiration Date	Please enquire.