

Product Datasheet

Mouse Serine Protease HTRA1 ELISA Kit (orb1210357)

Description	Mouse Serine Protease HTRA1, HTRA1 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	HTRA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,214
Uniprot ID	Q9R118
Solubility (25°C)	9600ng/L
Sensitivity	14.75ng/L
Entrez	56213
Expiration Date	Please enquire.