

Product Datasheet

Rat Mitogen Activated Protein Kinase 14 ELISA Kit (orb1210358)

Description	Rat Mitogen Activated Protein Kinase 14, MAPK14 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	MAPK14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,321
Uniprot ID	P70618
Solubility (25°C)	24ng/ml
Sensitivity	0.024ng/ml
Entrez	81649
Expiration Date	Please enquire.