

## Product Datasheet

### Human Coronin 2A ELISA Kit (orb1210369)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Human Coronin 2A, CORO2A<br>ELISA Kit  |
| <b>Reactivity</b>          | Human                                  |
| <b>Range</b>               | 0.05-15ng/ml                           |
| <b>Tested Applications</b> | ELISA                                  |
| <b>Target</b>              | CORO2A                                 |
| <b>Storage</b>             | 2-8°C                                  |
| <b>Note</b>                | For research use only                  |
| <b>Application notes</b>   | Sandwich                               |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates |
| <b>MW</b>                  | 59,763                                 |
| <b>Uniprot ID</b>          | <b>Q92828</b>                          |
| <b>Solubility (25°C)</b>   | 16ng/ml                                |
| <b>Sensitivity</b>         | 0.024ng/ml                             |
| <b>Entrez</b>              | <b>7464</b>                            |
| <b>Expiration Date</b>     | Please enquire.                        |