

Product Datasheet

Human Tyrosine Protein Phosphatase Non-receptor Type 11 ELISA Kit (orb1210370)

Description	Human Tyrosine Protein Phosphatase Non-receptor Type 11, PTPN11 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	PTPN11
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,011
Uniprot ID	Q06124
Solubility (25°C)	4000ng/L
Sensitivity	9.61ng/L
Entrez	5781
Expiration Date	Please enquire.