

Product Datasheet

Rat Regulating Synaptic Membrane Exocytosis Protein 1 ELISA Kit (orb1210372)

Description	Rat Regulating Synaptic Membrane Exocytosis Protein 1, RIMS1 ELISA Kit
Reactivity	Rat
Range	7-1500ng/L
Tested Applications	ELISA
Target	RIMS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	179,655
Uniprot ID	Q9JIR4
Solubility (25°C)	1600ng/L
Sensitivity	3.47ng/L
Entrez	84556
Expiration Date	Please enquire.