

Product Datasheet

Mouse Secretoglobin Family 1C Member 1 ELISA Kit (orb1210373)

Description	Mouse Secretoglobin Family 1C Member 1, SCGB1C1 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	SCGB1C1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,571
Uniprot ID	Q7M742
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	338417
Expiration Date	Please enquire.