

Product Datasheet

Human Cadherin EGF LAG Seven Pass G-type Receptor 2 ELISA Kit (orb1210374)

Description	Human Cadherin EGF LAG Seven Pass G-type Receptor 2, CELSR2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CELSR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	317,453
Uniprot ID	Q9HCU4
Solubility (25°C)	3200ng/L
Sensitivity	7.19ng/L
Entrez	1952
Expiration Date	Please enquire.