

Product Datasheet

Rat Ras-related Protein Rab-3A ELISA Kit (orb1210375)

Description	Rat Ras-related Protein Rab-3A, RAB3A ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	RAB3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,970
Uniprot ID	P63012
Solubility (25°C)	2400ng/L
Sensitivity	4.88ng/L
Entrez	25531
Expiration Date	Please enquire.