

Product Datasheet

Mouse Secreted Phosphoprotein 24 ELISA Kit (orb1210376)

Description	Mouse Secreted Phosphoprotein 24, SPP2 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/mL
Tested Applications	ELISA
Target	SPP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,136
Uniprot ID	Q8K1I3
Solubility (25°C)	1600ng/mL
Sensitivity	3.92ng/mL
Entrez	75396
Expiration Date	Please enquire.