

Product Datasheet

Mouse Scavenger Receptor Cysteine-rich Type 1 Protein M130 ELISA Kit (orb1210377)

Description	Mouse Scavenger Receptor Cysteine-rich Type 1 Protein M130, CD163 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	CD163
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	120,921
Uniprot ID	Q2VLH6
Solubility (25°C)	4800ng/L
Sensitivity	9.61ng/L
Entrez	93671
Expiration Date	Please enquire.