

Product Datasheet

Mouse R-Spondin-4 ELISA Kit (orb1210378)

Description	Mouse R-Spondin-4, RSPO4 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	RSPO4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,866
Uniprot ID	Q8BJ73
Solubility (25°C)	4800ng/L
Sensitivity	9.86ng/L
Entrez	228770
Expiration Date	Please enquire.