

Product Datasheet

Mouse R-Spondin-3 ELISA Kit (orb1210379)

Description	Mouse R-Spondin-3, RSPO3 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	RSPO3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,454
Uniprot ID	Q2TJ95
Solubility (25°C)	640ng/L
Sensitivity	0.97ng/L
Entrez	72780
Expiration Date	Please enquire.