

Product Datasheet

Human Homogentisate-1 ELISA Kit (orb1210380)

Description	Human Homogentisate-1, 2-dioxygenase, HGD ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	HGD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,964
Uniprot ID	Q93099
Solubility (25°C)	16ng/ml
Sensitivity	0.018ng/ml
Entrez	3081
Expiration Date	Please enquire.