

Product Datasheet

Rat Transcription Factor GATA-4 ELISA Kit (orb1210381)

Description	Rat Transcription Factor GATA-4, GATA-4 ELISA Kit
Reactivity	Rat
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	GATA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,597
Uniprot ID	P46152
Solubility (25°C)	48ng/ml
Sensitivity	0.044ng/ml
Entrez	54254
Expiration Date	Please enquire.