

Product Datasheet

Mouse Ribonuclease Pancreatic ELISA Kit (orb1210382)

Description	Mouse Ribonuclease Pancreatic , RNASE1 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	RNASE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,820
Uniprot ID	P00683
Solubility (25°C)	48ng/mL
Sensitivity	0.052ng/mL
Entrez	19752
Expiration Date	Please enquire.