

Product Datasheet

Mouse Rho-related GTP-binding Protein Rhoc ELISA Kit (orb1210383)

Description	Mouse Rho-related GTP-binding Protein Rhoc, RHOC ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	RHOC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,006
Uniprot ID	Q62159
Solubility (25°C)	2400ng/L
Sensitivity	5.24ng/L
Entrez	11853
Expiration Date	Please enquire.