

Product Datasheet

Mouse Rho-associated Protein Kinase 2 ELISA Kit (orb1210384)

Description	Mouse Rho-associated Protein Kinase 2, ROCK2 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	ROCK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	160,586
Uniprot ID	P70336
Solubility (25°C)	4800ng/L
Sensitivity	9.26ng/L
Entrez	19878
Expiration Date	Please enquire.