

Product Datasheet

Mouse Retinoic Acid Receptor Responder Protein 2 ELISA Kit (orb1210385)

Description	Mouse Retinoic Acid Receptor Responder Protein 2, RARRES2 ELISA Kit
Reactivity	Mouse
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	RARRES2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,350
Uniprot ID	Q9DD06
Solubility (25°C)	96ng/mL
Sensitivity	0.151ng/mL
Entrez	71660
Expiration Date	Please enquire.