

Product Datasheet

Mouse Resistin-like Beta ELISA Kit (orb1210386)

Description	Mouse Resistin-like Beta, RETNLB ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	RETNLB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,278
Uniprot ID	Q99P86
Solubility (25°C)	6400ng/L
Sensitivity	9.84ng/L
Entrez	57263
Expiration Date	Please enquire.