

Product Datasheet

Human Cartilage Associated Protein ELISA Kit (orb1210387)

Description	Human Cartilage Associated Protein, CRTAP ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CRTAP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,562
Uniprot ID	075718
Solubility (25°C)	4800ng/L
Sensitivity	9.17ng/L
Entrez	10491
Expiration Date	Please enquire.