

Product Datasheet

Rat L-lactate Dehydrogenase A Chain ELISA Kit (orb1210389)

Description	Rat L-lactate Dehydrogenase A Chain, LDH-A ELISA Kit
Reactivity	Rat
Range	0.05-30mU/ml
Tested Applications	ELISA
Target	LDH-A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,451
Uniprot ID	P04642
Solubility (25°C)	32mU/ml
Sensitivity	0.024mU/ml
Entrez	24533
Expiration Date	Please enquire.