

Product Datasheet

Mouse Regenerating Islet-derived Protein 3-gamma ELISA Kit (orb1210390)

Description	Mouse Regenerating Islet-derived Protein 3-gamma, REG3G ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	REG3G
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,307
Uniprot ID	009049
Solubility (25°C)	2400ng/L
Sensitivity	5.33ng/L
Entrez	19695
Expiration Date	Please enquire.