

Product Datasheet

Rat Aromatase ELISA Kit (orb1210391)

Description	Rat Aromatase, ARO ELISA Kit
Reactivity	Rat
Range	0.1-50ng/ml
Tested Applications	ELISA
Target	ARO
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,412
Uniprot ID	P22443
Solubility (25°C)	56ng/ml
Sensitivity	0.044ng/ml
Entrez	25147
Expiration Date	Please enquire.