

Product Datasheet

Mouse Receptor-interacting Serine Threonine-protein Kinase 3 ELISA Kit (orb1210392)

Description	Mouse Receptor-interacting Serine Threonine-protein Kinase 3, RIPK3 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	RIPK3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,322
Uniprot ID	Q9QZL0
Solubility (25°C)	6400ng/L
Sensitivity	9.05ng/L
Entrez	56532
Expiration Date	Please enquire.