

Product Datasheet

Human Clarin 3 ELISA Kit (orb1210393)

Description	Human Clarin 3, CLRN3 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CLRN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,321
Uniprot ID	Q8NCR9
Solubility (25°C)	4000ng/L
Sensitivity	9.34ng/L
Entrez	119467
Expiration Date	Please enquire.