

Product Datasheet

Mouse Receptor-interacting Serine Threonine-protein Kinase 1 ELISA Kit (orb1210394)

Description	Mouse Receptor-interacting Serine Threonine-protein Kinase 1, RIPK1 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/mL
Tested Applications	ELISA
Target	RIPK1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	74,854
Uniprot ID	Q60855
Solubility (25°C)	2400ng/mL
Sensitivity	4.55ng/mL
Entrez	19766
Expiration Date	Please enquire.