

Product Datasheet

Mouse Receptor Tyrosine-protein Kinase ERBB-3 ELISA Kit (orb1210395)

Description	Mouse Receptor Tyrosine-protein Kinase ERBB-3, ERBB3 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	ERBB3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	147,613
Uniprot ID	Q61526
Solubility (25°C)	48ng/mL
Sensitivity	0.054ng/mL
Entrez	13867
Expiration Date	Please enquire.