

Product Datasheet

Rat Keratin 19 ELISA Kit (orb1210396)

Description	Rat Keratin 19, KRT19 ELISA Kit
Reactivity	Rat
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	KRT19
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,636
Uniprot ID	Q63279
Solubility (25°C)	12.8ng/ml
Sensitivity	0.014ng/ml
Entrez	360626
Expiration Date	Please enquire.