

Product Datasheet

Mouse Ras-related Protein Rab-1A ELISA Kit (orb1210397)

Description	Mouse Ras-related Protein Rab-1A, RAB1A ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	RAB1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,678
Uniprot ID	P62821
Solubility (25°C)	3200ng/L
Sensitivity	7.23ng/L
Entrez	19324
Expiration Date	Please enquire.