

Product Datasheet

Rat Keratin 15 ELISA Kit (orb1210398)

Description	Rat Keratin 15, KRT15 ELISA Kit
Reactivity	Rat
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	KRT15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,870
Uniprot ID	Q6IFV3
Solubility (25°C)	16ng/ml
Sensitivity	0.017ng/ml
Entrez	287700
Expiration Date	Please enquire.