

Product Datasheet

Mouse Pyruvate Dehydrogenase Protein X Component ELISA Kit (orb1210399)

Description	Mouse Pyruvate Dehydrogenase Protein X Component, Mitochondrial, PDHX ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	PDHX
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,999
Uniprot ID	Q8BKZ9
Solubility (25°C)	3200ng/L
Sensitivity	7.33ng/L
Entrez	27402
Expiration Date	Please enquire.