

Product Datasheet

Mouse Ptb-containing ELISA Kit (orb1210400)

Description	Mouse Ptb-containing, Cubilin and LRP1-interacting Protein, PID1 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	PID1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,784
Uniprot ID	Q3UBG2
Solubility (25°C)	4800ng/L
Sensitivity	9.66ng/L
Entrez	98496
Expiration Date	Please enquire.