

Product Datasheet

Mouse Proto-oncogene Wnt-3 ELISA Kit (orb1210401)

Description	Mouse Proto-oncogene Wnt-3, WNT3 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	WNT3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,659
Uniprot ID	P17553
Solubility (25°C)	4800ng/L
Sensitivity	9.27ng/L
Entrez	22415
Expiration Date	Please enquire.