

Product Datasheet

Mouse Protocadherin Fat 3 ELISA Kit (orb1210402)

Description	Mouse Protocadherin Fat 3, FAT3 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	FAT3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	502,011
Uniprot ID	Q8BNA6
Solubility (25°C)	24ng/mL
Sensitivity	0.023ng/mL
Entrez	270120
Expiration Date	Please enquire.