

## Product Datasheet

# Rat Gamma-aminobutyric Acid Type B Receptor Subunit 1 ELISA Kit (orb1210403)

|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| <b>Description</b>         | Rat Gamma-aminobutyric Acid Type B Receptor Subunit 1, GABBR1 ELISA Kit |
| <b>Reactivity</b>          | Rat                                                                     |
| <b>Range</b>               | 10-2000ng/L                                                             |
| <b>Tested Applications</b> | ELISA                                                                   |
| <b>Target</b>              | GABBR1                                                                  |
| <b>Storage</b>             | 2-8°C                                                                   |
| <b>Note</b>                | For research use only                                                   |
| <b>Application notes</b>   | Sandwich                                                                |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                                  |
| <b>MW</b>                  | 111,534                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q9Z0U4</a>                                                  |
| <b>Solubility (25°C)</b>   | 2400ng/L                                                                |
| <b>Sensitivity</b>         | 4.99ng/L                                                                |
| <b>Entrez</b>              | <a href="#">81657</a>                                                   |
| <b>Expiration Date</b>     | Please enquire.                                                         |